FP12

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

```

move magpos 0.100 speed 1.00; wait ; section; idle 30;
section; move magpos 0.200 speed 1.00; wait ; section; idle 30;
section; move magpos 0.300 speed 1.00; wait ; section; idle 30;
section; move magpos 0.400 speed 1.00; wait ; section; idle 30;
section; move magpos 0.500 speed 1.00; wait ; section; idle 30;
section; move magpos 0.600 speed 1.00; wait ; section; idle 30;
section; move magpos 0.700 speed 1.00; wait ; section; idle 30;
section; move magpos 0.800 speed 1.00; wait ; section; idle 30;
section; move magpos 0.900 speed 1.00; wait ; section; idle 30;
section; move magpos 1.000 speed 1.00; wait ; section; idle 30;
section; move magpos 1.100 speed 1.00; wait ; section; idle 30;
section; move magpos 1.200 speed 1.00; wait ; section; idle 30;
section; move magpos 1.300 speed 1.00; wait ; section; idle 30;
section; move magpos 1.400 speed 1.00; wait ; section; idle 30;
section; move magpos 1.500 speed 1.00; wait ; section; idle 30;

```

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```

move magrot -80 speed 3.00; idle 6;
section;
move magrot -70 speed 3.00; idle 6;
section;
move magrot -60 speed 3.00; idle 6;
section;
move magrot -50 speed 3.00; idle 6;
section;
move magrot -40 speed 3.00; idle 6;
section;
move magrot -30 speed 3.00; idle 6;
section;
move magrot -20 speed 3.00; idle 6;
section;
move magrot -10 speed 3.00; idle 6;
section;
move magrot 0 speed 3.00; idle 6;
section;
move magrot 10 speed 3.00; idle 6;
section;
move magrot 20 speed 3.00; idle 6;
section;
move magrot 30 speed 3.00; idle 6;
section;
move magrot 40 speed 3.00; idle 6;
section;
move magrot 50 speed 3.00; idle 6;
section;
move magrot 60 speed 3.00; idle 6;

```

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FP12

Magnetic Tweezers Experiments to Probe DNA Mechanics

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Summary

Magnetic Tweezers

Forces and mechanical properties of biological macromolecules play an important role in biology, biophysics, and (bio)nanotechnology applications. Single-molecule magnetic tweezers are a biophysical method that can be used to investigate the mechanical properties of individual biomolecules, in particular DNA, RNA, or proteins. In the magnetic tweezers, molecules are fixed to a surface and bound at their other end to superparamagnetic beads approximately $1\text{ }\mu\text{m}$ in size. External magnetic field are used to exert defined stretching forces and torques on the beads and thus on the molecules. In this experiment, we investigate the behavior of DNA under different forces (so-called “force-extension curves”) and measure DNA supercoils, i.e., additional coils of the double-stranded DNA helix.

Contents

1	Introduction	1
2	Theory	4
2.1	Force-extension behaviour of molecules	4
2.1.1	Freely-jointed chain model (FJC)	4
2.1.2	Worm-like chain model (WLC)	6
2.2	Double-stranded DNA under force and torque	7
2.2.1	Force calibration	7
2.2.2	Applying twists to DNA molecules	9
3	Setup	11
3.1	Experimental setup	11
3.2	Software	13
3.3	Flow cell and attachment chemistry	13
4	Notes on this experiment	17
4.1	Questions for understanding	17
5	Measurements	19
5.1	Preparation of flow cell and instrument setup	19
5.2	Experiment 1: force-extension measurement	24
5.2.1	Force calibration from lateral fluctuations	24
5.3	Experiment 2: rotation-extension measurement	25
6	Analysis	26
6.1	Analysis for experiment 1	26
6.2	Analysis for experiment 2	26

1. Introduction

Mechanical properties and physical forces play important roles in many biological processes. An important example is DNA, which carries the genetic code and is increasingly used in a range of biotechnological applications. Our genetic information is stored in the sequence of DNA bases (A, C, T, and G). However, DNA is not just an abstract information carrier written in a four-letter alphabet, but also a physical polymer that is highly negatively charged and has distinct bending and torsional stiffnesses. These physical properties of DNA play crucial roles in the storage, readout and regulation of genetic information. For example, each of our cells contains more than 1 m of DNA, which must be stored in a $\approx 10\text{ }\mu\text{m}$ -sized nucleus. In addition, the binding of proteins to DNA molecules often induces a conformational change that can determine whether a gene is active or not. As all biological processes are influenced by the underlying physical properties of the molecules involved, the study of their mechanical properties is of great interest. There are several single-molecule techniques such as the atomic-force microscopy (AFM), optical tweezers (OT) or magnetic tweezers (MT), that can be used to measure mechanical parameters of molecules such as DNA or proteins at the single-molecule level.


Figure 1.1: (a) Principle of a MT setup: a superparamagnetic bead is attached to the DNA molecule using specific binding, while the other end is attached to the bottom of the flow cell through antibody binding. Permanent magnets are placed above the flow cell and exert a magnetic force on the bead. The stretching force and the rotation of the DNA molecule can be controlled by controlling the magnets position in translation and rotation. (b) Schematic illustrating the force range accessible to various single molecule techniques, extending up to several tens of picoNewtons (pN), a range relevant to probe biological processes at the molecular level. MT are particularly well suited for the precise measurement at very low forces.

In MT, a superparamagnetic bead is attached to the molecule of interest (here, DNA), while the other end is attached to the bottom of the flow cell. Two permanent magnets placed in a vertical configuration just above the flow cell exert a magnetic field that

induces a magnetic force on the magnetic bead and thus on the DNA molecule (**Figure 1.1**). The gradient of the magnetic field is nearly constant over typical molecule sizes, which implies that the magnetic force is constant over the field of view for a given distance of the magnets to the field cell. In contrast, OT and AFM measurements typically control molecular extension while the force fluctuates. Therefore, constant force measurements in OT and AFM require a feedback loop.

In MT, an inverted microscope and a CMOS sensor camera is used to track the magnetic bead's movement out-of-focus. From the bead's diffraction pattern, its three dimensional position (x , y , and z), and thus the DNA molecule's extension, can be determined. In addition, stretching force acting on the DNA molecule can be calibrated. Furthermore, DNA molecules can be twisted by rotating the external magnets.

In this advanced lab course, you will gain hands-on experience with the practical techniques of magnetic tweezers and the subsequent data analysis. The experiments focus on investigating the mechanical properties of DNA. You will learn the principles of a conventional magnetic tweezers setup, perform standard force-extension measurements of double-stranded DNA, and conduct rotational experiments to observe DNA supercoiling.

Learning objectives

After completing the experiment, you will be able to . . .

1. explain the operating principle of magnetic tweezers and distinguish them from optical tweezers,
2. determine the force acting on a magnetic bead from its Brownian fluctuations,
3. describe force-extension data of a single DNA molecule using the FJC and WLC models,
4. determine the persistence length and Kuhn length from experimental data,
5. explain how torsional stress and supercoiling arise in DNA,
6. determine a buckling point from a rotation-extension curve,
7. estimate the torsional stiffness of DNA from the measurement data,
8. discuss experimental deviations from the idealized models.

2. Theory

Double-stranded DNA molecules are typically much longer than their diameter of ≈ 2 nm. The length of double-stranded DNA depends on the number of base pairs, with each base pair contributing ≈ 0.34 nm. For the 3 billion base pairs of the human genome, the total contour length is ≈ 1 m. In this work, we use a 21,000 bp DNA molecule, which is approximately $7 \mu\text{m}$ long—more than 1,000-fold longer than its width. For such chain-like molecules, statistical models can be applied to describe their response to stretching forces.

In this laboratory course, you will perform force-extension measurements on double-stranded DNA, which means to apply various stretching forces and observe the DNA's tether extension. In a second experiment you will apply turns to the DNA molecule and study its behavior upon analyzing changes in the DNA's length. This chapter covers briefly the theoretical background that is needed to do this laboratory course.

Table 2.0.1: Notation of the key quantities used in this chapter.

Symbol	Meaning	Unit
L_c	DNA contour length	μm
z	measured extension	μm
F	stretching force	pN
L_p	persistence length	nm
l_K	Kuhn length	nm
C	torsional stiffness	pN nm^2
n	number of applied turns	1

2.1 Force-extension behaviour of molecules

In order to describe the response of constant forces acting on biological macromolecules, we would like to introduce two typical theoretical models, the freely-jointed chain model (FJC) and the worm-like chain model (WLC). The first model (FJC) considers the molecule as a chain, consisting of rigid segments, while the other model (WLC) considers the molecule as an isotropic rod that is continuously flexible.

2.1.1 Freely-jointed chain model (FJC)

The FJC model considers a polymer of total length L that consists of N monomers each with the same length l (Kuhn-length):

$$L = N \cdot l$$

$\vec{R}$ is the end- to- end vector of the molecule, describing the "actual" size of the molecule in its environment. All segments are described by a vector $l \cdot \vec{r}_i$. The FJC model assumes that there are no interactions between the monomers and the bond between two vectors $\vec{r}_i$ and $\vec{r}_{i+1}$ is perfectly flexible (and the chain, therefore, performs a random walk), see **Figure 2.1**. All conformations of the polymer are equally likely to occur in solution. If an external force $\vec{F} = f \cdot \vec{e}_z$ acts on the molecule, the potential energy E_{pot}^α of the molecule in a particular angular conformation α is:

$$E_{pot}^\alpha = - \sum_{i=1}^N \vec{F} \cdot l \cdot \vec{r}_i = -l \cdot f \sum_{i=1}^N \cos \theta_i \quad (2.1.1)$$

Where θ is the angle between one monomer and f the force acting in z- direction.


Figure 2.1: In the FJC model a polymer is considered to consist of N rigid segments with length l . Also shown in this figure is the definition of the end- to- end vector $\vec{R}$, which is spanned from the first segment to the last segment of the polymer.

The response (extension) of a molecule to a force acting on it, can be derived through statistical mechanics. A canonical ensemble is a statistical ensemble to represent all possible states of a mechanical system that is in thermal equilibrium. In a canonical ensemble the volume V , the number of particles N and the temperature T are considered constant. The canonical partition sum Z can be defined using the following formula and inserting equation (2.1.1) leads to the following result:

$$Z = \sum_{\alpha} e^{\frac{-E_{pot}^\alpha}{kT}} = \sum_{\alpha} \prod_{i=1}^N e^{\frac{fl \cos \theta_i}{kT}} = \left(\sum_{\alpha} e^{\frac{fl \cos \theta}{kT}} \right)^N \quad (2.1.2)$$

Where k is the Boltzmann constant and T is the absolute temperature.

Due to the fact that the angle is continuous, the sum can be replaced by an integral, integrating over the whole space:

$$Z = \left(\int d\Omega \cdot e^{\frac{fl \cos \theta}{kT}} \right)^N = \left(\frac{2\pi kT}{fl} \sinh\left(\frac{fl}{kT}\right) \right)^N \quad (2.1.3)$$

Inserting this to the definition of the *Gibbs free energy* $G = -kT \cdot \ln(Z)$ the (mean) extension z of the molecule due to the force $\vec{F}$ acting on it can be calculated:

$$\langle z \rangle = -\frac{\delta G}{\delta f} = L \cdot \left(\coth\left(\frac{fl}{kT}\right) - \frac{kT}{fl} \right) \quad (2.1.4)$$

2.1.2 Worm-like chain model (WLC)

In comparison to the FJC model, the WLC model considers the molecule as a continuously flexible chain, see **Figure 2.2**.


Figure 2.2: The worm-like chain is comparable to the classical mechanical bending of a rod, see black box on the top left. The molecule of length L_c is parametrized by points s on the molecule. Vector $\vec{r}(s)$ points to the chain whereas $\vec{t}(s)$ is the tangential vector to point s .

The chain consists of segments s with length ds , with $s \in [0, L_c]$, where L_c is the contour length of the molecule. $\vec{t}(s)$ are the tangent vectors to each point s . If a force $\vec{F}$ is acting on the molecule, it acts on each tangential vector. Therefor the formula of the energy E^{WLC} ends up in an integral, integrating over the whole length of the molecule:

$$E = \int_0^{L_c} ds \cdot \vec{t}(s) \cdot \vec{F} = f \cdot \int_0^{L_c} ds \cdot \vec{t}(s) \cdot \vec{e}_z \quad (2.1.5)$$

Generally the worm-like chain is solved in analogy to the freely-jointed chain: first, one defines the energy for one conformation (E_α^{WLC}) of the polymer and second, one calculates the partition sum. In order to integrate over all conformations, one has to perform a path integral $D(\vec{t}(s))$, integrating over all tangent vectors in point s . Here, we keep the mathematics to a minimum. Consider that the polymer behaves like a bar that is bend by a force $\vec{F}$ (classical mechanics), see **Figure 2.2**, black box. The energy to deform a bar is proportional to $\frac{1}{R^2}$, where R is the radius of the bending circle. Taking into account that $(\frac{\delta \vec{t}(s)}{\delta s})^2 = \frac{1}{R^2}$ we get:

$$E_\alpha^{WLC} = \int_0^{L_c} ds \cdot \left(\frac{A}{2} \cdot \left(\frac{\delta \vec{t}(s)}{\delta s} \right)^2 \right) \quad (2.1.6)$$

Where $A = kT \cdot L_P$. L_P is called persistence length and defines the rigidity of the chain. Segments, smaller than the persistence length behave like a rigid rod. Segments

larger than the persistence length are flexible and can be bend. An exact solution of the WLC model is hard to derive, as one has to perform the path integral:

$$\int_{\vec{t}(0)}^{\vec{t}(L)} D(\vec{t}(s)) \cdot \exp\left(-\frac{L_P}{2}\right) \cdot \int_0^L ds \left(\frac{\delta \vec{t}(s)}{\delta s}\right)^2 \quad (2.1.7)$$

Instead of an exact numerical solution one often uses an approximate interpolation formula, which is sufficiently accurate to analyse experimental data:

$$\frac{f L_P}{kT} = \frac{z}{L} + \frac{1}{4 \cdot (1 - \frac{z}{L})^2} - \frac{1}{4} \quad (2.1.8)$$

Where f is the force acting on the molecule, L_P is the persistence length, z is the extension and L_C is the total length of the molecule.

A key point here is that it requires energy (supplied by the external -magnetic, in this case- force) to stretch a polymer in the FJC or WLC model, despite the fact that the model assumes that there are no interaction between the segments. The reason is that there are many conformation of the polymer where the ends are close together and only few (in the extreme case of being fully stretched only one) configuration with the ends far apart. Stretching the polymer is, therefore, entropically unfavourable and this situation is sometimes called an “entropic spring”.

2.2 Double-stranded DNA under force and torque

In this chapter, we explain how the stretching force on the bead is calibrated in a magnetic tweezers setup and how DNA reacts to induced torques.

2.2.1 Force calibration

In a magnetic tweezers setup, magnets exert a magnetic field and thus a magnetic force on a magnetic bead. This bead is attached to one end of the DNA molecule (via multiple biotin-streptavidin linkages), while the other end is attached to the bottom of the flow cell via antibody binding, see **Figure 2.3**. The bead transfers the force to the molecule: it can be stretched and twisted simply by controlling the magnets position in translation and rotation.

The stretching force is calibrated using the fluctuations of the magnetic bead in solution, which are caused by Brownian motion. The setup can be seen as an inverted pendulum, see **Figure 2.3**.

The force ($\vec{F}_m$) on the bead is pointing in z-direction. Due to Brownian motion the bead fluctuates around its origin. A restoring force ($\vec{F}_R$) forces the bead always back to its initial position. Using the angular deviation Θ and the bead’s transverse fluctuations δx the force is defined:

$$\vec{F}_R = \sin \Theta \cdot \vec{F}_m \approx \frac{\delta x}{\langle z \rangle} \cdot \vec{F}_m \quad (2.2.9)$$


Figure 2.3: The magnets exert a force $\vec{F}_m$ on the magnetic bead. The force can be calibrated using the deviation δx of the bead from its original position, which is due to Brownian motion. The length of the molecule is labeled with L , the angular deviation is called Θ and $\vec{F}_R$ is the restoring force acting on the bead.

We find that the bead is in a harmonic potential with spring constant $\kappa = \frac{F_m}{\langle z \rangle}$. By the equipartition theorem, each degree of freedom has $\frac{1}{2}k_B T$ energy in thermal equilibrium:

$$\frac{1}{2} \cdot \kappa \cdot \langle \delta x^2 \rangle = \frac{1}{2} k_B T \quad (2.2.10)$$

Inserting the spring constant and solving the equation for F_m , we can calculate the force acting on the bead:

$$F_m = \frac{k_B T \langle z \rangle}{\langle \delta x^2 \rangle} \quad (2.2.11)$$

where, k_B is the Boltzmann constant, T is the absolute temperature, $\langle z \rangle$ is the mean extension of the DNA tether, and $\langle \delta x^2 \rangle$ is the variance of the bead's lateral fluctuations.

In other words, the force can be calibrated by measuring the the x-y movement of the bead for a sufficiently long time. To ensure reliable force-extension measurements, the experimental protocol must account for the system's characteristic relaxation time (τ) to determine the total measurement time (t) required for statistical accuracy. The relaxation time defines time scale of the bead undergoing thermal fluctuations within the harmonic trap. It is calculated as the ratio of the friction coefficient (γ) to the trap's spring constant (k):

$$\tau = \frac{\gamma}{k}$$

The trap stiffness k in the magnetic tweezers is given by the ratio of the applied force to the measured extension of the DNA tether $\frac{F}{\langle z \rangle}$. Consequently, higher forces result in a shorter relaxation times and thus faster bead fluctuations.

To reach a desired statistical accuracy (ϵ), the measurement time must be significantly

longer than the relaxation time (τ), long enough to sample these fluctuations sufficiently. You can think about the accuracy improving with $\sqrt{N}$, where N is the number of independent measurements. In this case, each measurement of duration τ can be thought of as another independent measurement. For forces above 1 pN, the minimum measurement time τ_{measure} is estimated as:

$$\tau_{\text{measure}} \approx \frac{12\pi^2\eta R l_0}{F\epsilon^2}$$

where, η is the fluid viscosity, R is the bead radius, and l_0 is the DNA contour length.

This relationship implies that data collection time must be increased when working with longer DNA tethers or simply when aiming for higher precision (smaller ϵ , typically 5% or 10%) to ensure statistically meaningful force–extension measurements.

2.2.2 Applying twists to DNA molecules

In addition to applying stretching forces, MT can also twist DNA molecules. In a typical measurement, the DNA is held at a constant stretching force and is either overwound or underwound by n turns. DNA is a right handed molecule and behaves (for forces larger 1 pN) different upon the direction of applied turns. Positive turns twist DNA with its helical nature, while negative turns twist the molecule against its intrinsic helicity. The response of DNA to applied turns depends on the stretching force. For forces smaller 1 pN, DNA behaves symmetric according to the applied number of turns, see **Figure 2.4A** (left column). For a small number of turns the molecule's extension stays constant, while the molecular torque increases linearly. Beyond the buckling transition (indicated with blue squares for 0.9 and 1.1 pN) the molecule's extension decreases linearly, while the molecular torque stays constant. At the buckling transition it is energetically more favorable to form DNA plectonemes as indicated by coarse-grained simulations in 2.4B in contrast to staying overwound. In the post-buckling regime every turn of the molecule leads to the formation of a DNA supercoil. For higher stretching forces (> 1 pN) the respond of DNA upon the direction of applied turns is asymmetric. While in general the buckling point shifts to a larger number of applied turns with higher stretching forces, there is no buckling point in negative direction, see **Figure 2.4A** (right column). Instead torque-induced melting occurs at -10 pN·nm. Torsional stress is hereby released by local basepair breaks.

At the buckling point (n_b) the DNA's extension is reduced due to the formation of plectonemic supercoils. The change in the end-to-end distance corresponds to the part of DNA that forms DNA supercoils. A simple model for DNA buckling and supercoiling is as follows. In the pre-buckling regime, when the DNA remains approximately straight, the torque and torsional energy increase as the molecule is twisted more and more with the number of applied turns:

$$E_{\text{torsion}} = 2\pi n \cdot \Gamma \quad (2.2.12)$$

$$\Gamma = \frac{C}{L} \cdot 2\pi n \quad (2.2.13)$$

where Γ is the torque and C is the torsional stiffness.


Figure 2.4: (A) Extension (top) and torque vs. applied number of turns (bottom) for different stretching forces. Buckling points (n_b) are indicated for 0.9 pN and 1.1 pN in the extension vs. turns plot with a blue square. Torque measurements were obtained using magnetic torque tweezers, a specialized type of magnetic tweezers. (B) Coarse-grained DNA simulations show torsionally relaxed DNA, underwound DNA (upon torque-induced melting at -10 pN·nm) and overwound DNA with the formation of DNA plectonemes.

At the buckling point (n_b) the twist energy $E_{torsion}$ of the molecule corresponds to the energy that is needed to form such a loop (E_{loop}). E_{loop} has two contributions: the work W that is needed to shorten the DNA and E_{bend} to bend the DNA strand.

$$E_{torsion}(n_b) = E_{loop} = W + E_{bend} \quad (2.2.14)$$

In our simple model, one supercoil is approximated as a circular loop with radius r . The work, needed to reduce the length of the DNA molecule against the force F is

$$W = -2\pi r F \quad (2.2.15)$$

while the energy for the bending to perform one loop is (see the WLC section)

$$E_{bend} = \int_0^{2\pi r} ds \frac{A}{2} \frac{1}{r^2} = \pi A \frac{1}{r} \quad (2.2.16)$$

Minimizing the looping energy with respect to r , one can derive the optimal radius r^* of the loop

$$r^* = \sqrt{\frac{A}{2F}} \quad (2.2.17)$$

Having the radius, one can calculate the buckling torque Γ_b and the number of turns n_b where the buckling takes place.

$$\Gamma_b = \sqrt{2AF} = \sqrt{2L_P k T F} \quad (2.2.18)$$

$$n_b = \frac{L}{2\pi C} \cdot \sqrt{2AF} \quad (2.2.19)$$

3. Setup

This chapter describes the experimental setup, the software used to control the instrument, and the fabrication of the flow cells used in this experiment.

3.1 Experimental setup

In MT, the molecule of interest (here, DNA) is attached to the surface of a flow cell and on the other side to a superparamagnetic bead (**Figure 3.1**). A pair of permanent magnets placed above the flow cell exert magnetic fields, which in turn induce magnetic forces on the superparamagnetic bead. The stretching force can be controlled by changing the permanent magnets position, i.e. translating it closer and further away from the flow cell surface. The molecules' rotation can be controlled by rotating the magnets. The flow cell is mounted on a flow cell holder, such that the flow cell is placed on top of an oil immersion objective. The objective can be moved up and down very precisely using a piezo, which is in particular needed to perform a "look up table" (LUT) for z-tracking. To control liquid exchange, a peristaltic pump is connected to the flow cell. A LED monochromatically illuminates the flow cell through the gap (1 mm) of the magnet assembly. Images of the illuminated area are taken by a CMOS camera and transferred to a computer where image analysis is performed. The DNA tether is not directly observed in the MT, instead the bead positions are tracked and used for data analysis.

Measurements of beads are performed out-of-focus. The diffraction pattern of the beads is used in order to track the tether extension along z (**Figure 3.2**). A calibration of the bead's diffraction pattern is performed prior to the measurement. For the calibration, the magnets are brought close to the flow cell to exert large stretching forces which reduce the fluctuations of the bead itself. The position of the objective can be controlled by a piezo. In order to create a LUT, the objective's position is moved in sub- μm steps and the bead's diffraction pattern is captured. To get precise information in vertical direction, the software interpolates between each picture and creates the LUT. During the measurement, the piezo is kept at a fixed position and the current diffraction pattern of the beads are compared to the LUT, which returns the z-position of the bead and, therefore, the extension of the molecule. For the lateral position of the bead, the center of the diffraction rings are calculated via cross-correlation.

To reduce the effect of mechanical drift, which can disturb the measurements in lateral as well as in vertical direction, reference beads are tracked. These reference beads are polystyrene beads that are stuck to the flow cell surface. To correct for drift, the recorded trace of a reference bead is subtracted from the DNA-bound magnetic beads under study.


Figure 3.1: Overview of a MT setup. A pair of permanent magnets is placed on top of the flow cell. The light of the LED illuminates the sample from above. The light is directed to a camera via a mirror and a lens.


Figure 3.2: Tracked pictures of magnetic beads. The bead's diffraction pattern changes depending on its distance from the objective's focal plane. This feature is used to determine the extension of bead-bound DNA.

You will be asked to perform a force calibration of the MT setup as well as force-extension measurements for DNA molecules. **Figure 3.3** illustrates the forces as a function of the height of the magnets above the flow cell used in this laboratory course and can be used to write measurement protocols.


Figure 3.3: Example of a force vs. magnet position calibration for the MT setup used in this lab course with MyOne beads. The solid line is a double-exponential fit to the measured values. Inset: the data on a semi-log plot with logarithmic y-axis, and linear x-axis

3.2 Software

The software to control the instrument is written in LabVIEW. The piezo stage and both motors responsible for magnet translation and rotation, can be controlled using the software. The peristaltic pump to exchange fluids has to be used manually.

The software consists of three separated windows (see **Figure 3.4** and **3.5**). A live image of the field of view of the flow cell is shown in “BeadTrackerMain_modV1.vi” window. Basic setup configurations, selection of beads and the LUT are made in this window. Motors and the piezo stage are controlled in “MotorUI.vi”. “ExperimentProgramUI_modV1.vi” contains all options to start and control the measurements. The bead positions are tracked in real time and saved in text files including the motor positions.

3.3 Flow cell and attachment chemistry

The flow cell (**Figure 3.6a**) provides the measurement chamber for the magnetic tweezers. Two microscope cover slips separated by a soft melted parafilm spacer provide a $\approx 25 \mu\text{L}$ microfluidic chamber. The top cover slip has two holes (inlet and outlet) for fluid exchange in the flow cell. The bottom coverslip is functionalized with (3-glycidyloxypropyl)trimethoxysilane followed by incubation with anti-digoxigenin for DNA attachment. The linear 21 kbp DNA has two “handles”, i.e. DNA segments at the end that serve as attachment points. Each handle is approximately 600 bp long and contains several nucleotides with biotin (biotin handle) or digoxigenin groups (digoxigenin handle). The biotin handle binds to a streptavidin coated $1 \mu\text{m}$ diameter superparamagnetic bead and the digoxigenin labeled end of the DNA binds to the anti-digoxigenin coated surface of the flow cell (**Figure 3.6b**). Often, $1 \mu\text{m}$ diameter polystyrene beads are used as reference beads. Note that torque and twist measurement on DNA are only possible if the double-stranded DNA does not have any “defects” called nicks and if the attachment is via multiple biotin and streptavidin moieties at the ends. To avoid nonspecific binding




Figure 3.4: “BeadTrackerMain_modV1.vi” is the main window of the software. It displays a real-time image of the field of view and provides option to select beads and perform a LUT.

of the bead to the surface, the cover slip is passivated with bovine serum albumin (BSA).


(a) The “ExperimentProgramUI.modV1.vi” window allows users to define the folder path, create protocols, and start measurements.

(b) With “MotorUI.vi”, the piezo-stage as well as the two motors for the magnets can be controlled manually.

Figure 3.5: “ExperimentProgramUI.modV1.vi” and “MotorUI.vi” to start measurements and to control motor positions, respectively.

(a) **Schematic illustration of a flow cell.** Two glass coverslips separated by a parafilm spacer form a microfluidic chamber for MT measurements. Flow cell assembly is placed in a flow cell holder and is mounted in the optical stage of the MT setup. The flow cell is connected to a peristaltic pump for controlled liquid exchange.

(b) **Schematic illustration of a bead-tether-surface construct.** Anti-digoxigenin is bound to the bottom coverslip by epoxysilane chemistry. The paramagnetic bead has a streptavidin coated surface. A DNA construct containing digoxigenin labels on one end and biotin on the other binds to the surface and the bead, respectively. For passivation, the bottom coverslip is coated with BSA.

Figure 3.6: Flow cell construction and attachment chemistry for typical MT measurements.

4. Notes on this experiment

Please submit your files using the filename structure **FP12_GRXXX_Jahr_Monat_Tag**. All important information can be found at <https://www.uni-augsburg.de/de/fakultaet/mntf/physik/groups/exp4/teaching/fp/>.

Preparation for the practical day

To ensure a smooth workflow and successful experimentation in the lab, we strongly ask you to **arrive well prepared for the experiment day**. A thorough engagement with the theoretical foundations (especially polymer models such as FJC and WLC as well as DNA supercoiling) and with the individual steps of the measurement protocol is mandatory. To support this and to promote general understanding, you will find a preparatory question catalogue at the end of this chapter.

Use of digital devices in the S1 laboratory

Please note that the practical lab is located in a safety level 1 (S1) area. For reasons of work safety and contamination prevention, we ask you to **avoid using digital end devices** (e.g. iPads, tablets, laptops) to read the script in the lab whenever possible.

- **Personal notes:** If you would like to use handwritten notes or your own annotations, please print the script in advance and bring it in paper form.
- **Printouts on site:** If you do not wish to print your own copy, a printed version of the instructions or step-by-step protocol is available directly in the lab at the chair.

4.1 Questions for understanding

Use the following questions to check your preparation for the pre-test and the measurements:

1. Operating principle & measurement setup:

- Why are superparamagnetic rather than ferromagnetic beads used in magnetic tweezers?
- Via which biochemical bonds is the DNA molecule attached to the surface of the flow cell and to the magnetic bead?

- How is the Z-position (extension) of the beads determined optically even though the DNA molecule itself is not visible under the microscope? What is the purpose of creating a look-up table (LUT)?

2. Force calibration & polymer models:

- How can the tensile force acting on the molecule be calibrated using Brownian motion and the equipartition theorem?
- How do the basic assumptions of the *Freely-Jointed-Chain* (FJC) and *Worm-Like-Chain* (WLC) models differ?
- What is meant by the persistence length L_P and the Kuhn length l ? Why is the stretching of DNA in the low-force regime described as an „entropic spring“?

3. DNA supercoiling & topology:

- What is meant by the buckling point and what happens mechanically to the DNA beyond this threshold?
- Why does double-stranded DNA respond asymmetrically to positive and negative turns when $F > 1$ pN?
- What is a DNA *nick* and how does a nicked DNA become noticeable in the rotation experiment?

4. Experimental procedure:

- What are reference beads used for during the measurement?
- Which three preliminary tests are carried out before the actual data acquisition to ensure that only single, intact, and torsionally constrained DNA tethers are measured?

Use of aids and AI tools

As in any scientific work, the principle of good scientific practice also applies to the protocol and evaluation for this practical course: independent work and the tools used must always be clearly and transparently identified. If you use generative AI tools (such as ChatGPT or similar language models) for text creation, formulation support, or code generation, this use must be explicitly stated in the report and explained at which points and to what extent the respective tool was used. Purely independent work must remain clearly recognizable to the assessors.

5. Measurements

This chapter provides you step-by-step guidance for the experiments. Please read the experimental protocol carefully as outlined here.

5.1 Preparation of flow cell and instrument setup

The flow cell has already been pre-coated with anti-digoxigenin and stored overnight. To attach the DNA, prepare a DNA-bead solution and inject it into the flow cell. After incubation, flush the microfluidic chamber with buffer (we use PBS or “phosphate buffered saline”, a frequently used buffer that mimics the pH and salt concentration found in cells) to remove unbound beads, and finally insert the magnets. The magnet offset step has also been performed beforehand.

Note: The magnet offset defines the reference “zero” position used for the experiments. All subsequent magnet positions are measured relative to this offset. To determine the effective offset, start with the magnet position slightly above the flow cell and then decrease it carefully in step with small increments (0.01 mm) and carefully monitor the beads’ movement. The magnet position at which bead motion due to the magnets touching the top surface of the flow cell is first detected is recorded as the magnet offset. This defines the magnet’s zero position relative to the flow cell.

A step-by step protocol is provided below. Before starting with the following protocol, remove the magnet holder and the motor has to be set to “Magnet all up” (≈ 18.9 mm).

- Flush the flow cell with 600-800 μL 1x PBS using the peristaltic pump.
- Flush the flow cell with 300 μL of BSA in PBS (100 μL at a time).
- Fill the inlet with PBS and incubate for atleast 1.5 hours.
- Flush flow cell with 600-800 μL 1x PBS (200 μL at a time) and wrap the inlet of the flow cell holder with parafilm.
- Cleaning the beads: add 2 μL of well mixed bead solution (MyOne superparamagnetic 1 μm beads) to 198 μL 1x PBS and mix by brief vortexing in a LoBind microfuge tube.
- Place the tube in a magnetic holder to draw the beads to one side, wait 20 seconds, and then remove the liquid.
- Redo the washing step two times.

- Add 50 μL 1x PBS.

Note: Briefly vortex the solution to make sure all the beads are in the buffer, before adding DNA.

- Add 2 μL DNA construct to the bead solution and incubate a couple of minutes with the microfuge tube kept in horizontal position. Avoid mixing with a pipette (results in DNA nicks).
- Fill up with 48 μL 1x PBS to have overall volume of 100 μL .
- Fill the chamber with ca. 50 μL DNA-beads and wait for 20 minutes.

Note: If no DNA tether is observed at the end, use the remaining 50 μL of the DNA-bead solution to restart from this step.

- Rise the chamber with 1x PBS until the number of freely moving beads is greatly reduced.
- Stop the pump and insert the magnets. Apply high force (≈ 5 pN, that is 0.6 mm magnet distance) and rotations (e.g. alternative rotation between 0 and 50 turns with 5 turns/s) to loosen the nonspecifically bound beads.
- Remove the magnet and repeat rinsing step with PBS until no unbound beads remain.
- Mount the magnets back in the holder.

Note: From this point, the magnets should never have a larger distance than 10 mm from the flow cell for more than 20 seconds to avoid sticking of the beads to the microfluidic chamber surface.

In the next step, select the beads that are bound to a single, coilable DNA molecule.

- Create a folder in the computer to save the location of the beadlist. Follow this path: ExperimentProgramUI_modV1.vi $\rightarrow$ Path and Settings $\rightarrow$ ‘Choose/auswählen’.
- Move the piezo-stage up/down (window “MotorUI.vi”) to get the beads in focus (ideally one should find two foci, one for beads attached to DNA and one for beads unspecifically attached to the surface of the flow cell).
- For a coarse check of DNA-bound beads, bring the magnets little closer to the flow cell (e.g. 2 mm) and look for beads that show real-time lateral fluctuations (due to the presence of magnetic field) within the field of view in the BeadTrackerMain_modV1.vi window. Alternatively, the DNA-bound beads can be focused using the piezo, as their patterns will differ from those of the beads on the surface.
- Further, search for three beads which are not moving, i.e. unspecifically attached to the surface, and whose diffraction pattern are not overlapping. Remember their positions. These beads can be used as reference beads.

Note. This is only valid when the reference beads are selected from the stuck myOne beads. While using polystyrene beads, reference beads are typically chosen before adding the DNA.

- Change the ROI to 120.
- Press “select beads” and click “remove all”. Press “add” and click first on the three beads chosen before as reference beads. Now, click on every other DNA-bound beads whose diffraction patterns are not overlapping. Make sure that the squares that the program puts around the beads is fully visible in the field of view. Do not use beads that are too close to the edge of the field of view.
- Press “close” to close the window and “save” the beadlist.

For z-tracking, a LUT has to be calibrated for every single bead:

- In the “Tracking” panel, set the number of steps to 200 with a step size of $0.1\ \mu\text{m}$. This defines the total Z range i.e. $20\ \mu\text{m}$.
- Move magnets to a distance of 0.6 mm to apply a high force ($\approx 5\ \text{pN}$), which is necessary to avoid large fluctuations of the beads during the LUT-measurement.
- In the “MotorUI.vi” window, adjust the focal length to define the lower limit of the LUT. This should be done by playing with the piezo (objective) distance with respect to the beads and monitoring the diffraction pattern originating from the out-of-focus beads.
- Click on “Tracking” in ‘BeadTrackerMain_modV1.vi’ and set current piezo Z as the start piezo Z (lower limit of LUT).
- Press “Build ZLUT” to start the LUT measurement and wait until the software has finished. Save the ZLUT.

Note. In case the software crashes, one can load the bead list and the ZLUT again.

The setup is now ready to track beads. To ensure that tracked beads are attached to a single and coilable DNA tether, we perform three tests prior to the actual measurements:

- Choose a folder to save the measurements and a name for the measurement in window “ExperimentProgramUI_modV1.vi”.
- The 1st test is to check DNA tethers that have expected contour length. For that, the magnet distance to the flow cell is varied from 0.6 mm to 10 mm, which corresponds to $\approx 5\ \text{pN}$ and close to zero force for the MyOne beads used in this experiment.

Note: At a magnet height of 0.6 mm, the DNA tethers are stretched to $> 95\%$ of their contour length whereas at 10 mm (close to zero force) the beads drop the bottom surface, corresponding to zero tether extension.

- Start the measurement with the button “Run experiment” in the window “ExperimentProgramUI_modV1.vi”.
- A new window will pop up. Wait until the camera reacts (frame number in this window) before clicking on any other button.
- Move the motor to 10 mm distance, wait 10 seconds and drive back to 0.6 mm.

- Stop measurement (“abort experiment”).
- Activate “Trace & LUT” in window “BeadTrackerMain_modV1.vi” and click through the traces. Only beads with tether lengths close to the expected contour length of 21 kbp DNA ($\approx 7 \mu\text{m}$) are selected for the further measurements. (In this screening, do not consider the beads number 0, 1 and 2 which are reference beads).

Note: Subtract a reference bead from the traces of DNA-bound beads to see their z-extension more clearly.

- Change the measurement name for the 2nd test.

Note: The 2nd measurement is to check if any bead is bound to multiple DNA tethers and to select those beads which are tethered by a single double stranded DNA molecule. This is achieved by the asymmetric mechanical response of DNA upon underwinding and overwinding at different stretching forces.

- Check that the motor is at a distance of 0.6 mm ($\approx 5 \text{ pN}$).
- Start the measurement.
- Rotate the motors to -70 turns (at 2 turns/s) i.e. negative turns under high tension.
- Rotate back to 0 turns at the same stretching force ($\approx 5 \text{ pN}$).
- Stop the measurement.
- Check the traces for beads who do not behave as expected and delete them with “Discard bead”.

Note: Bead showing a decrease in tether length is assumed to be bound by multiple DNA molecules.

- Change the measurement name for the 3rd test.

Physical background: A DNA nick is a break in only one of the two DNA strands. It removes the torsional coupling of the DNA molecule: applied turns can relax at the break and are no longer stored as torsional stress throughout the molecule. In the rotation experiment, a nick therefore produces a markedly different rotation curve, in particular the absence or strong reduction of the characteristic tether shortening after buckling.

Note: This is to check whether any torsionally unconstrained tether (with breaks or nicks) is present among the selected beads.

Note: The 3rd test is typically performed together with the 2nd test so that both tests can be plotted together. This allows for simultaneous evaluation of the two characteristics and simplifies creating an overview of the beads when recording the results.

- Rotate + 70 turns at a distance of 3.7 mm ($\approx 0.5 \text{ pN}$) and rotate back to 0 turns keeping the same magnet distance (i.e. same stretching force). Check the traces and discard beads which do not behave as expected.

Note: A bead attached by a single, coilable DNA molecule shows buckling and a corresponding decrease in the tether length, in contrast to nicked DNA molecules.


- Only beads those pass all these tests will be considered for the experiments.

5.2 Experiment 1: force-extension measurement

5.2.1 Force calibration from lateral fluctuations

For each magnet distance used in the protocol, record not only the z-position but also the bead's lateral x-y fluctuations. Determine the force from the variance of the x-fluctuation using the equipartition theorem. Plot the resulting force values against magnet distance. Fit the force-magnet-position curve with a double exponential function and use it to determine the force values for the force-extension and rotation-extension measurements.

Track the beads' position (x, y, z) by varying the magnet position with respect to the flow cell. An example of a measurement protocol is shown in **Figure 5.1**. Start with the highest force in the protocol. Before running the script, press “Simulate script” in the “ExperimentProgramUI.vi”.


```

move magpos 0.100 speed 1.00; wait ; section; idle 30;

section; move magpos 0.200 speed 1.00; wait ; section; idle 30;
section; move magpos 0.300 speed 1.00; wait ; section; idle 30;
section; move magpos 0.400 speed 1.00; wait ; section; idle 30;
section; move magpos 0.500 speed 1.00; wait ; section; idle 30;
section; move magpos 0.600 speed 1.00; wait ; section; idle 30;
section; move magpos 0.700 speed 1.00; wait ; section; idle 30;
section; move magpos 0.800 speed 1.00; wait ; section; idle 30;
section; move magpos 0.900 speed 1.00; wait ; section; idle 30;
section; move magpos 1.000 speed 1.00; wait ; section; idle 30;
section; move magpos 1.100 speed 1.00; wait ; section; idle 30;
section; move magpos 1.200 speed 1.00; wait ; section; idle 30;
section; move magpos 1.300 speed 1.00; wait ; section; idle 30;
section; move magpos 1.400 speed 1.00; wait ; section; idle 30;
section; move magpos 1.500 speed 1.00; wait ; section; idle 30;

```

Figure 5.1: An example of a script for the force-extension measurement. Move the position of the magnet in steps, the value following “magpos” is the absolute distance of the magnets from the flow cell. The number following “idle” defines the time to record frames in seconds. The speed value should be always used as given in the example (z-motor = 1.0).

5.3 Experiment 2: rotation-extension measurement

Here, test the DNA tether extension upon turning the magnets. Therefore, measure rotation curves from - 70 to + 70 turns in steps of 2 turns at forces of 5.0 pN, 2.0 pN, 1.5 pN, 1.0 pN, 0.7 pN and 0.5 pN. Start with the highest force. It is recommended to start a new measurement for every force for easier analysis. A measurement time of 6 s for each turn is sufficient. Alternatively, run the protocol shown in **Figure 5.2** at different stretching forces.

The force values used in this experiment are obtained from the force–magnet–position calibration determined in Experiment 1.

```

move magrot -80 speed 3.00; idle 6;
section;
move magrot -70 speed 3.00; idle 6;
section;
move magrot -60 speed 3.00; idle 6;
section;
move magrot -50 speed 3.00; idle 6;
section;
move magrot -40 speed 3.00; idle 6;
section;
move magrot -30 speed 3.00; idle 6;
section;
move magrot -20 speed 3.00; idle 6;
section;
move magrot -10 speed 3.00; idle 6;
section;
move magrot 0 speed 3.00; idle 6;
section;
move magrot 10 speed 3.00; idle 6;
section;
move magrot 20 speed 3.00; idle 6;
section;
move magrot 30 speed 3.00; idle 6;
section;
move magrot 40 speed 3.00; idle 6;
section;
move magrot 50 speed 3.00; idle 6;
section;
move magrot 60 speed 3.00; idle 6;
section;

```

Figure 5.2: An example of a script for the rotation-extension measurement. The value following “magrot” defines the absolute values of the rotation motor. The number following “idle” defines the time to record frames in seconds. The speed values should be always used as given in the example (rotation motor = 3.0). Every “section” will save the frame numbers in a text file, which makes it easy to find motor movements for analysis.

6. Analysis

6.1 Analysis for experiment 1

- Determine the force for each magnet distance from the variance of the lateral x-fluctuation using the equipartition theorem. Generate a plot of force against magnet distance, similar to **Figure 3.3**, and fit a double exponential function to your data. Use this calibration curve to determine the force values for the force-extension and rotation-extension measurements.
- Plot the force-extension curve of the DNA molecule and fit both the FJC model and the WLC model. Determine the Kuhn length from the FJC fit and the persistence length from the WLC fit. How many base pairs does the persistence length determined from the WLC fit correspond to? What is the physical significance of this length? Compare your values to literature. Compare the quality of the FJC and WLC descriptions in the investigated force range. In which ranges do the models show systematic deviations from the data? Discuss possible physical causes.

6.2 Analysis for experiment 2

- Plot the extension of the DNA molecule as a function of the number of turns for each applied force.
- Explain the differences observed in the plots for different forces.
- Determine the approximate size of a single supercoil.
- Determine the buckling point n_b and the buckling torque Γ_b for all forces and determine the torsional stiffness of DNA, comparing your results with Strick *et al.* [4].

Bibliography

- [1] S. Cocco, J. F. Marko, and R. Monasson. Theoretical Models for Single-Molecule DNA and RNA Experiments: From Elasticity to Unzipping. *Comptes Rendus. Physique*, 3(5):569–584, 2002.
- [2] D. Dulin. An introduction to magnetic tweezers. In I. Heller, D. Dulin, and E. J. G. Peterman, editors, *Single Molecule Analysis: Methods and Protocols*, volume 2694 of *Methods in Molecular Biology*, pages 375–401. Humana Press, New York, NY, 3 edition, 2023.
- [3] J. Lipfert, J. W. J. Kerssemakers, T. Jager, and N. H. Dekker. Magnetic Torque Tweezers: Measuring Torsional Stiffness in DNA and RecA-DNA Filaments. *Nature Methods*, 7(12):977–980, 2010.
- [4] T. R. Strick, J. F. Allemand, D. Bensimon, A. Bensimon, and V. Croquette. The Elasticity of a Single Supercoiled DNA Molecule. *Science*, 1996.
- [5] I. D. Vilfan, J. Lipfert, D. A. Koster, S. G. Lemay, and N. H. Dekker. Magnetic tweezers for single-molecule experiments. In P. Hinterdorfer and A. van Oijen, editors, *Handbook of Single-Molecule Biophysics*, pages 371–395. Springer US, New York, NY, 2009.